PolyBranch HCR™ RNA FISH Kit(Animal)


Catalog Number:HKR12D

【Product Information】

Product Name	Catalog No	20T	Storage	Shelf Life
Wash Buffer (powder, dissolve in 10L ddH ₂ O)	HKR12-1	10L	Short-term:4℃ Long-term:-20℃	12 months
Proteinase K (100x)	HKR12-2	20uL		
HRP-Mouse Anti-Digoxin (100x)	HKR12-3	20uL		
TSA Chromogenic Solution 570nm (Optional)	HKR12-4	2mL		
Blocking Buffer	HKR12-5	2mL		
Probe:	HKR12-6	2mL		
Amplification Probe	HKR12-7	2mL		

【Product Description】

RNA hybridization kit, suitable for in situ hybridization of RNA larger than 300bp, with good specificity and high sensitivity. The reagent kit contains multiple pairs of combination probes. When two probes bind to adjacent positions of the target sequence, a small segment of the top sequence can trigger a branched HCR to form a large molecular nucleic acid aggregate. Each nucleic acid strand is labeled with digoxin, and the antibody recognizes digoxin for color development.

【Probe Information】

【Tissue Fixation】

Tissue Type	Animal	Plant	Frozen Samples	Cell Climb Slides	Cells
Treatment	Fix at RT for 12h, paraffin embed	Vacuum fix for 1h, RT fix for 12h, paraffin embed.	Dehydrate in 15% sucrose at 4℃ for 8h, then in 30% sucrose at 4℃ for 8h, OCT embed.	Fix at 4℃ for 2h.	scrape off cells, fix in 4% PFA at 4℃ for 2h, wash with PBS, agarose embed.
Type	mRNA	lncRNA	circRNA	miRNA	rRNA
Treatment	Fix at RT for 12h	Fix at RT for 12h (<300bp: 24h)	Fix at RT for 12h	Fix at RT for 12h	Fix at RT for 12h (<300bp: 24h)

【Storage and Shipping】

Ship on wet ice; store at -20℃ for long-term or at 4℃ for short-term use. Shelf life: 6

months.

【Protocol】

1. Deparaffinization and Rehydration

Immerse slides sequentially in: Xylene I (15 min) → Xylene II (15 min) → Xylene III (15 min) → 100% ethanol (10 min) → 90% ethanol (10 min) → 80% ethanol (10 min) → 70% ethanol (10 min) → Rinse with distilled water.

2. Enzyme Repair

After slides are completely dry, draw a hydrophobic circle around the tissue using a histology pen (recommended: HKR14P In Situ Hybridization Pen). Place slides horizontally in a hybridization oven or humidified chamber. Add 100 μ L of Proteinase K repair solution (1X) onto the tissue and incubate at 37°C for 30 min (15 min for cell samples). Rinse with distilled water to stop the reaction.

3. Blocking

Remove excess liquid from slides. Add 100 μ L of Blocking Buffer per slide and incubate at 37°C for 30 min. Wash once for 5 min. **Wash steps: Place slide racks in a wash tank, add wash buffer (ensure samples are submerged), and shake at 60 rpm for 5 min.**

4. Probe Hybridization

Remove excess liquid from slides. Add 100 μ L of probe per slide and incubate at 37°C for 3 h or overnight (maintain humidity to prevent drying). Wash 5 times, 5 min each. **Wash steps: As described above.**

5. Probe Amplification

Remove excess liquid from slides. Add 100 μ L of Amplification Probe per slide and incubate at 37°C for 1.5 h (maintain humidity). Wash 5 times, 5 min each. **Wash steps: As described above.**

6. HRP-Mouse Anti-Digoxin

Remove excess liquid from slides. Add 100 μL of HRP-Mouse Anti-Digoxin (1X) per slide and incubate at 37 °C in a humidified chamber for 40 min. Wash 5 times, 5 min each.

Wash steps: As described above.

7. Chromogenic Reaction

Remove excess liquid from slides. Add 100 μL of TSA Chromogenic Solution per slide and incubate at RT for 10 min. Rinse with distilled water to stop the reaction.


8. DAPI Staining

Add 50 μL of DAPI staining solution per slide, incubate in the dark for 5 min, rinse with distilled water, and mount with anti-fade mounting medium.

【Precautions】

1. For research use only.
2. Wear lab coats and disposable gloves for safety.

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